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COMUNICADO

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RELEASE

TO THE ATTENTION OF THE DISTRIBUTORS AND USERS OF THE HOISTS
WINNMOTION 175, BASICA E 130 AND HOP E 150

With the motivation of increasing user safety, on June 9, 2016, the assembly process of the Winnmotion 175, Basica E130 and HOP E150 hoists was improved.

This January of 2018 we have gone back to improve some details of the process of assembly, quality control, information of use and training on the central point of suspension of the hanger (safety nut). Everything properly collected in the instructions for use.

In these instructions we also include the ANNUAL INSPECTION OF PATIENT HOISTS where we indicate the control points of the hoist to ensure its perfect condition.

This annual maintenance (ISO 10535, Annex B) is mandatory and the best way to avoid incidents such as those that have already occurred on 2 occasions in France due to negligence in the maintenance of the devices named Winnmotion 175.

We appreciate your trust and respect for these measures necessary for the proper use and better experience and user satisfaction.

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Cecilio CASAS
Responsable Qualité et SAV
Winnicare Spain
ccasas@winnicare.es

Attachment: Updated instructions for use.

WINNCARE SPAIN SLU
B-96221718
Carretera Masía del Juez 37b
46909 - Torrent, Valencia (Spain)
Tfno. (+34) 96 156 55 21
FAX : +34 9 61 56 03 83
EMAIL : sat@winncare.es
<http://www.winncare.es>

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To must read before use, to return to the user and keep

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
UK

The lock nut and spring pin system should in no way be any attempt at loosening or adjustment. This system, mounted at the factory, allowing securing the scourge. Disassembly or intentional release of this set spring pin and lock nut can lead to loosening of the scourge that may cause the patient to fall.

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