
Rehosoft®

≤ 120 kg
≤ 265 lb

10 cm

MD

WINNCARE FRANCE, Miniac
ZONE ACTIPOLE
2 RUE DE SAINT-COULBAN
35540 MINIAC-MORVAN - FRANCE

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Compensation for disability.

Indication

Going to the toilet in people with locomotor and neuromuscular disorders (decreased force in the knees, hips, ankles; altered postural balance).

Contraindication

Patient with a maximum weight greater than that tolerated by the raised toilet seat.

Target patient and user group

Dependent patients, amputees, patients who underwent hip surgery, patients with difficulty sitting and/or standing up.


Indicate undesirable side effects

Any serious incident related to the device should be notified to the manufacturer and the competent authority of the Member State where the user and/or the patient is established.

Inform the competent authority if you consider or have reason to believe that the device presents a serious risk or that it is a fake product.

2. COMPOSITION OF THE MEDICAL DEVICE

- ❶ PVC shell containing injected polyurethane foam.
- ❷ Toilet lid.

3. CLINICAL BENEFIT, PERFORMANCE

Device performance characteristics

Allows the disabled or older person to sit and stand up, reducing the risk of falls.

Expected clinical benefit

Compensation for a decrease in muscle strength, joint aging, and multisensory organization of postural control.

4. INSTRUCTIONS FOR USE

Device installation

The product is ready to be installed. Simply place it on the toilet bowl.

Maintenance

- > Between each patient.
- > Processes: bio-cleaning, steam cleaning.
- > Products: detergents and surface disinfectants.
- > Compatible with bleach solution of up to 1000 ppm.

5. WARNINGS AND PRECAUTIONS USE

Precautions for use

- > When you lower the toilet lid, move your fingers away from the raised toilet seat to prevent them from jamming.
- > Before each use, check that the product is not damaged and that it is properly assembled.
- > Make sure that the raised toilet seat is placed in the middle of the bowl.

Warnings

- > Incorrect use of the raised toilet seat may result in injury or damage.
- > Do not modify the device.
- > Never stand on the raised toilet seat.
- > Do not sit on the toilet lid.
- > A "sit-to-sit" transfer from the wheelchair to the raised toilet seat should be assisted.
- > In order to avoid the risk of a patient fall, before performing a "stand-to-sit" transfer, the caregiver needs to verify that the device is centred and stable.

Lifetime information

The device may be reused within the initial performance maintenance limits during its life cycle. Its estimated lifespan is 4 years.

Accessories / spare parts

Only the toilet seat lid and spare parts that are supplied by the manufacturer can ensure the safety and performance of the complete product.

6. DISPOSAL

Do not dispose of the product in nature outside of the designated areas. Follow the recycling channels that are in place in your country.

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Siotec AG, Sägestrasse 75, CH-3098 Köniz
Mangar International Limited (WINNOCARE group)
Cubro Limited. 149 Taurikura Drive, Tauriko,