

DAGONE ROOM



LIGNES CONTEMPORAINES



FURNITURE FOR THE DAGONE ROOM



BOARDS

- DIMENSIONS:
Head: W1005 x H561 mm
Foot: W1005 x H641 mm
- ▶ Curved hand rail
 - ▶ Inside boards are in melamine-faced MDF



DESK

- DIMENSIONS:
W1085 x D548 x H758 mm



BEDROOM TABLE

- DIMENSIONS:
W885 x D548 x H758 mm
- ▶ 1 drawer
 - ▶ Not compatible with Bridge



BEDSIDE CABINET

- DIMENSIONS:
W455 x D447 x H681 mm
W548 x D560 x H814 mm with refrigerator
- ▶ 2 rear mobility castors
 - ▶ 2 adjustable feet at the front
 - ▶ 1 drawer with removable tray
 - ▶ Hinge door hung to the left
 - ▶ 1 shelf with bottle space



DRESSER

- DIMENSIONS:
W885 x D548 x H758 mm
- ▶ 3 drawers



TABLE WITH DRESSER

- DIMENSIONS:
W1426 x D548 x H758 mm
- ▶ 3 right-hand drawers



SINGLE-DOOR WARDROBE

- DIMENSIONS:
W685 x D548 x H1793 mm
- ▶ left-hand hinge
 - ▶ 1 shelf
 - ▶ 1 clothes rail



TWO-DOOR WARDROBE

- DIMENSIONS:
W1135 x D548 x H1793 mm
- ▶ 1 full-width top shelf
 - ▶ 3 half-shelves
 - ▶ 1 right-hand clothes rail

DAGONE	FURNITURE AND BOARD COLORS		SURROUND, EDGING AND HAND RAIL COLORS	MATERIAL
	STANDARD COLORS	OTHER COLORS		
BOARDS	Inside finish in Light Beech. Ral 1001 Lacquered hand rails and edging	Choice of melamine finishes from our color chart. Veneered edging	Edging and surround in ral 1001 as standard.	MDF hand rail and melamine faced MDF inside
FURNITURE				Unit in melamine-faced 19 mm chipboard / 1 mm ABS edging Tops in MDF with 19 mm melamine-face / Lacquered edging

1-year warranty*. *Under normal conditions of use and excluding wearing parts - please see instructions of use

* Please note: since these ranges are in MDF, you are strongly advised against choosing the Light Beech varnish; you may opt for the Light Beech lacquer instead. (this is because light varnish does not go well on MDF parts).

Photos not contractually binding

DAGONE Atmosphere
Furniture and boards: Dagone, Taupe, ral 9016
Bridge: PVC Atoll, structure ral 9016

FP000016 / 01.2018



WINNCARE Group - 4, Le Pas du Château - 85670 Saint-Paul-Mont-Penit - France
Tél : +33 (0)2 51 98 55 64 - Fax : +33 (0)2 51 98 59 07 - E-mail : info@medicatlantique.fr - www.winnicare.fr
In the interest of constantly improving its products, the WINNCARE Group reserves the rights to modify the technical specifications of any product without prior notice
Class I medical device in accordance with European Directive 93/42/EEC

